






MECCA HEALTHCARE PRIVATE LIMITED



CATALOGUE

OUR MISSION

Our mission is to be an organization whose employees strive for growth and values the contribution made by its founding members, through always keeping in mind the health and safety of the users of its products and thus satisfying the growing needs of the healthcare industry and its customers.

OUR VISION

To be an organization whose people & products exude CARE, COMPASSION & TRUST towards its customers, business partners and the society at large.

A BRIEF NOTE ABOUT MECCA INDUSTRIES

"Mecca Group of Industries", is a company set up in the year 1977 catering to the healthcare industry. Mecca Industries is among the first Indian company to begin mfg. Intravenous Infusion Equipments.

Our organization has three state of the art mfg. plants located at Kalol Dist; Gandhinagar & Chhatral Dist. Mehsana in Gujarat and Boranada Dist; Jodhpur, in Rajasthan.

Today we are an CE MARKED PRODUCTS, ISO 13485:2016 and WHO GMP compliant concern supplying to varied sectors of the healthcare industry in India. Recently Mecca has started developing 510K approval US FDA facility at its Kalol plant and intends to procure the certification by 2020.

Our products on an annual rate contract basis are supplied to State Governments of Gujarat, Madhya Pradesh, Maharashtra, Bihar and Tamil Nadu in India.

Mecca group is also exporting its CE marked products to all the CIS countries, Middle Eastern Nations, South East Asian countries, African Nations and Latin American countries since 1991 through contract manufacturing and merchant export.

We are also working on a loan license contract mfg. basis with Indian arms of two of the largest German companies of healthcare sector in the world, **M/s. B. Braun India Pvt. Ltd. and M/s. Fresenius Kabi India Pvt. Ltd. & M/s Hollister Inc. USA. We are working on Make In India Initiative to supply specialised devices for advanced drug delivery therapy with M/s Hetero Labs and M/s Cipla.**

Besides we are holding loan licence of Claris Otsuka Ltd (Formerly Core Parenterals associated since 1988). In past our co. has worked with almost all the Intravenous Infusion mfg. cos in India prominent being Sarabhai Industries, Cadila Pharmaceuticals, Tide Pharmaceuticals (a unit of Torrent Pharma.), Wockhardt Pharma, Intas Pharmaceuticals, Cipla Limited, Albert David Ltd. et al.

To know more about our organization and its line of work we suggest you to visit. www.mhplindia.com

Mecca Healthcare (P) Limited, INDIA

Plant Address;

UNIT - 1, 99 - 102, GIDC, KALOL (N.G.), Dist; Gandhinagar, Gujarat

UNIT - 2, L-1202, GIDC Chattral, Dist.; Mehsana, Gujarat

UNIT - 3, F-252, III Phase, Boranada, Jodhpur, Rajasthan

www.mhplindia.com

AN ISO 13485 : 2016 & ISO 9001 : 2015 CERTIFIED COMPANY

WHO GMP COMPLIANT MFG. PLANT

CE MARKED PRODUCTS

Infusion Products

Needle Free Connector



1. Swab the top of the Needle Free Connector with an alcohol pad for at least three seconds or per your facility's protocol.



2. Insert the syringe or connect I.V. tubing directly to the Needle Free Connector and turn clockwise.



3. Turn the syringe or I.V. set counterclockwise to remove. Replace the Needle Free Connector according to your facility's protocol.

Features: & Specifications

- 1.Transparent housing material: Polycarbonate or Copolyester.
- 2.Metal-free & MRI-compatible.
- 3.Latex-free.
- 4.ISO 10993 compliant.
- 5.Multiple insertion at least 100 times/day.
- 6.Priming Volumn: 0.09ml.
- 7.Perfect fl owing rate: 350 ml/min under one meter water pressure, tested by Mecca Technical Team.
- 8.Produced under by GMP: Mecca is an ISO9001, ISO13485 EU MDRA Govt. of India CDSCO registered manufacturing facility.

Infusion Products

Y-port Valves



Material: Polycarbonate or Copolyester
Size: ID 4.01mm, 3.6mm

Tube End Valves



Material: Polycarbonate or Copolyester
Size: ID 2.0mm, 2.5mm, 2.8mm, 3.7mm, 4.1mm

WE SPECIALIZE IN

MECCA INFUSION DEVICES

Lightproof Photosensitive IV Administration Infusion Set (Micro Infusion Set)

Product Specifications

- This Infusion set is prepared from Lightproof / Photosensitive polymer materials. Angioplast provides this Photosensitive sets available in Yellow and Brown color and are obtainable with or without Precision Flow Controller [Dial Type] 5 ml to 250 ml/hrs.



Features of our Lightproof Photosensitive Infusion Set:-

- Latex Free.
- Soft, Clear and Flexible Tubing.
- Sterile, Non-Toxic and Pyrogen-Free.
- Sterilized by E.O. (Ethylene Oxide)



- 88" IV Administration Set Vented, Vented Spike 20 Drops-ml with Guard/Brown Color Drip Chamber 60mm/Fluid Filter 15 Micron/ Flow Regulator/Precision Flow Controller 5-250ml-hr. /Latex Free 'Y'-Site/Rotating Luer Lock Male Adaptor with Cap/3.0 x 4.1mm PVC Tubing Brown Colour Length 30+140+30cm/Sterile/Paper Poly Pack/ Quantities per Carton = 300 pcs



- 88" IV Administration Set Vented, Vented Spike 20 Drops-ml with Guard/Yellow Color Drip Chamber 60mm/Fluid Filter 15 Micron/ Flow Regulator/Latex Free Flash Bulb/Rotating Luer Lock Male Adaptor with Cap/ 3.0 x 4.1mm Yellow Color PVC Tubing Length 30+140+30cm/Sterile/Paper Poly Pack/ Quantities per Carton = 300

DCS

INNOVATIVE "AIR STOP" FEATURE KEY BENEFITS

1. Avoids "Air Infusion" when IV Bottle gets empty & maintains a constant fluid level inside Infusion line so no air can get inside the patient body.
2. Unique Technology "Precision-Filter" acts as a barrier for particulate matter 15 microns & below.
3. Offers "No-Hassle", "Risk-Free" & "Smooth Switch-Over" to next IV Bottle.

INNOVATIVE "PRIME STOP" FEATURE KEY BENEFITS

1. "Auto Fill" of Infusion Set with "No Dripping" of IV Fluid on hands, floors & bedding at time of IV Infusion.
2. Offers ease in Multiple Solution preparations and saves time.
3. Maintains a closed system before the commencement of IV Therapy to a patient.

- * Fitted with Safety Roller Controller offering secure disposal of IV Set after use.
- * Fitted with Rotating Luer Lock for easy fitment to an IV Catheter.
- * Fitted with a Slip Clamp offering ease in temporary stop-over of IV Therapy.
- * Fitted with a Y-Site for Extra Medication

KEY BENEFITS

	Mecca IV Set	Normal IV Set
Air Embolism Risk Free	Yes	No
Time and Energy Savings in Multiple Solution Preparation	Yes	No
Prevention of Infusion Line getting Dry	Yes	No
Smooth Switch Over to next IV Bottle	Yes	No
Maintenance of a closed system before IV Therapy	Yes	No

A Technical Collaboration with ACME UK Inc.
Limited, London, Great Britain

Manufactured At: 99, GIDC, Kalol, Gandhinagar,
Gujarat, India, Pin: 382721

Contact for a Demo: +91 999 869 9090





- 84" IV Administration Set Vented, Vented Spike 20 drops-ml with Guard/Drip Chamber 70mm/Fluid Filter 15 micron/Flow Regulator/ Dial Flow Regulator/ Latex Free 'Y'-Site/0.2 Micron Filter/Rotating Luer Lock Male Adaptor with Cap/ 3.0mm x4.1mm PVC Tubing-Length 60+50+30+15cm/EO Sterile/Paper Poly Pack / Quantities per Carton = 300 pcs



- it-01a1 82" IV Administration Set Non-Vented, Non-Vented Spike 10 drops-ml with Guard/Drip Chamber 70mm/Flow Regulator/Latex Free 'Y' Site/0.2 Micron Filter/Rotating Luer Lock Male Adaptor with Cap/3.0mm x4.1mm PVC Tubing-Length 115+30+15cm/EO Sterile/Paper Poly Pack / Quantities per Carton = 300 pcs



- 82" IV Administration Set Vented, Vented Spike 20 drops-ml with Guard/Drip Chamber 70mm/Fluid Filter - -15 micron/Flow Regulator/Latex Free 'Y'- Site/0.2 Micron Filter/Rotating Luer Lock Male Adaptor with Cap/ 3.0mm x 4.1mm PVC Tubing-Length 115+30+15cm/EO Sterile/Paper Poly Pack / Quantities per Carton = 300 pcs



- 12" Standard Bore Extension Set, Female Luer with Cap/0.2 Micron Filter/Slide Clamp/Latex Free 'Y'- Site/Rotating Luer Lock Male Adaptor with Cap/3.0mm x4.1mm PVC Tubing-Length 10+10+10cm/EO Sterile/ Paper Poly Pack / Quantities per Carton = 400 pcs



- 12" Standard Bore Extension Set, Female Luer with Cap/1.2 Micron Filter/Slide Clamp/Latex Free 'Y' Site/Rotating Luer Lock Male Adaptor with Cap/3.0mm x4.1mm PVC Tubing-Length 10+10+10cm/EO Sterile/ Paper Poly Pack / Quantities per Carton = 400 pcs

Dial Flow IV Infusion Set (Micro Infusion Set)

Product Specifications

- Dial flow IV Infusion set comprises of a device that can be "dialed" to a specific infusion rate. They are ideally suited for gravity infusions and available in both vented (for bottles and bags) and non-vented (for Plastic bottles).



Some characteristics of our Dial flow IV Infusion set are:-

- Sharp spike for better penetration
- Reliable flow rate accuracy across the full adjustable scale Latex Free.
- Soft Clear Flexible Tubing.
- ASterile, Non-Toxic and Pyrogen Free.
- Rotating Lock to simplify connection to patient.
- Pinch clamp to facilitate immediate shut- off when necessary.
- Custom configurations available upon request.



- Vented Spike with Guard/20 drops-ml/Drip Chamber 60mm/Fluid Filter 15 micron/Flow Regulator/ Dial Flow Regulator/ Latex Free 'Y'-Site /Rotating Luer Lock Male Adaptor/ 2.7mm x 3.9mm PVC Tubing Length 50+130+15cm/EO Sterile/Paper Poly Pack / Quantities per Carton = 240 pcs



- Vented Spike with Guard/60 drops-ml/Drip Chamber 60mm/Fluid Filter 15 micron/Micro Drip Tube/Flow Regulator/Dial Flow Regulator/ Latex Free 'Y'- Site /Rotating Luer Lock Male Adaptor/2.7mm x3.9mm PVC Tubing-Length 50+130+15cm/EO Sterile/Paper Poly Pack / Quantities per Carton = 240 pcs

Infusion Set With BACK CHECK VALVE


Product Specifications

- The device is efficient in pumping any liquid into our body. The tubing of the device comes in variable lengths to accommodate various needs.



Some salient features of our Infusion set for Infusion Pump are as follows:-

- Available in vented (for bottles and bags) and non-vented (for plastic bottles)
- Sharp spike for better penetration.
- ISO 8536-4 compliant spike ISO 594/1 & ISO 594/2 compliant end-connection.
- Latex Free.
- Soft, Clear and flexible Tubing.
- Sterile, Non-Toxic and Pyrogen Free.
- Sterile Pack - ETO (Ethylene Oxide) Sterilization validate according to Harmonized Standards.



- 62" IV Administration Set Non-Vented, Non-Vented Spike 10 Drops-ml with Guard/Drip Chamber 60mm/Fluid Filter-15 micron/ Latex Free 'Y'Site/Flow Regulator/Luer Tip Adaptor with Cap/3.0mm x 4.1mm PVC Tubing-Length 135+15cm/EO Sterile/Paper Poly Pack / Quantities per Carton = 300 pcs



- 62" IV Administration Set Vented, Vented Spike 20 Drops-ml with Guard/Drip Chamber 70mm/Fluid Filter-15 micron/ Latex Free 'Y'-Site/Flow Regulator/Luer Lock Adaptor with Cap/3.0mm x 4.1mm PVC Tubing-Length 135+15cm/EO Sterile/Paper Poly Pack / Quantities per Carton = 300 pcs

Paclitaxel Infusion Set

(Micro Infusion Set)

Product Specifications

- Self and fast priming of adult I.V. set to curtail nursing time. Automatic air-venting in any position loss of flow due to air locks. Able to withstand a minimum of 45psi (3.1 bar) pressure. ISO 8536-4 compliant Spike & ISO 594/1 and 594/2 compliant- end - connection. Bacterial retention for over 96 hrs. Flow. Options available with PVC/DEHP - Free. Options available with Filter media 1.2 / 0.2 micron. EO Sterile, Non - Toxic and Pyrogen Free.



- 76" Infusion for Paclitaxel Administration, Vented Spike 20 drops-ml with Guard/PVC-DEHP- Free Drip Chamber 60mm/ Fluid Filter 15micron/Flow Regulator/0.2 Micron Filter/Rotating Luer Lock Male Adaptor with Cap/3.0mm x4.1mm PVC-DEHP Free Tubing-Length 140+15cm/EO Sterile/Paper Poly Pack Quantities per Carton = 300

CONTACT US

Corporate Office

Mecca Industries, B-605, Ratnakar Nine Square,
Opposite ITC Narmada, Judges Bungalow Road,
Near Mansi Crossroads, Vastrapur,
Ahmedabad, Gujarat, India, Pincode: 380015
Tel; +91 79 4800 5989

UNIT 1 ADDRESS

102, Mecca Industries Compound,
G.I.D.C.Kalol (N.G.), Gandhinagar,
Gujarat, India,
Pincode: 382725
Tel; +91 2764 221020

UNIT 2 ADDRESS

F-252, III Phase, Boranada, Jodhpur,
Rajasthan, India, Pincode: 342012
Tel; +91 8320107430

Branch Office

L-1202, GIDC Chattral, Dist.
Mehsana, Gujarat India

REGD. OFFICE

408, Spectrum, Near G.P.O., Relief
Road, Ahmedabad-01, Gujarat
India

OUR UK REPRESENTATIVE'S ADDRESS

ACME UK INC LTD

25 Leeming Road, Borehamwood,
Hertfordshire WD64EB
United Kingdom

Website: www.acmevision.com